

Automated evaluation of disease on ^{68}Ga -DOTATATE PET/CT images for long-term lesion tracking on ^{177}Lu -DOTATATE therapy


Timothy G Perk, Amy J Weisman, Ojaswita Lokre, Brayden Schott,
Victor Santoro-Fernandes, Robert Jeraj, Steve Y Cho, Scott B Perlman




Patient Population

- Retrospective study of 25 patients
- Sequential ^{68}Ga -DOTATATE PET/CT imaging during ^{177}Lu -DOTATATE therapy
- 81 images examined
- 12 patients had 3 or more images (range 2-7)
- Images analyzed using TRAQinform IQ[®] technology (AIQ Solutions)


Example ^{68}Ga -DOTATATE PET/CT Image




Early Heterogeneity of Response (^{177}Lu -DOTATATE)



14	New
2	Increasing
16	Stable
19	Decreasing
9	Disappeared



9	New
10	Increasing
7	Stable
2	Decreasing
11	Disappeared



- 24/25 patients exhibited heterogeneous response of ROI at the first imaging timepoint
- 11/25 patients have over half of their disease increasing or new

Summary

- AIQ Technology can automatically identify and track ROI on ^{68}Ga -DOTATATE PET/CT images over the course of ^{177}Lu -DOTATATE treatment
- All 25 patients showed heterogeneous ROI response on at least one timepoint
- 58% of patients (7/12) showed an escalation in the number of ROI that were decreasing or disappeared over time
- Understanding the response of each ROI can help inform treatment strategies